

## ESS Lock N' Load Handles and Syringes

The easiest way to collect undisturbed, measured soil cores for field preservation techniques for EPA Method 5035.

- Handle locks at 5 and 10 gram settings with one easy turn.
- Beveled edge syringe is stronger than cut syringes and fits in the neck of 40 mL glass vials.
- Turn handle and push to dispense soil into the vial. Soil dispensing is done without having to remove the syringe.
- Strong, Accurate, Easy and Economical

Lock N' Load Handle and Syringes can also be used in the lab. Great for subsampling from brass sleeves or soil jars to 40mL vials for auto-sampler extraction.



**TTT Environmental**  
Instruments and Supplies

*The preferred source for instrument  
Rentals, Sales, Service and Supplies!*

4201 B Street  
Anchorage, AK 99503  
(907) 770-9041

[www.tttenviro.com](http://www.tttenviro.com)

25420 74<sup>th</sup> Avenue S.  
Kent, WA 98032  
(253) 373-9041

[info@tttenviro.com](mailto:info@tttenviro.com)