# 2100P and 2100P IS Portable Turbidimeter

## Features and Benefits

#### **Laboratory Quality in a Portable Unit**

The Hach 2100P and 2100P IS Portable Turbidimeters offer a level of performance previously possible only with laboratory instruments. Microprocessor-controlled operation and Hach's unique Ratio  $^{\text{TM}}$  optics bring great accuracy, sensitivity, and reliability to field and in-plant testing.

#### Two Models for Specific Requirements

- 2100P Turbidimeter—Get fast, accurate turbidity testing in the field or the lab, over a wide range of samples.
   Compliant with USEPA Method 180.1 design criteria.
- 2100P IS Turbidimeter—Designed to meet international standards that mandate measurement using an LED light source.

#### **Two-detector Optical System**

The two-detector optical system compensates for color in the sample, light fluctuation, and stray light, enabling analysts to achieve laboratory-grade performance on a wide range of samples, even under difficult, onsite conditions.



The Hach 2100P and 2100P IS Portable Turbidimeters bring laboratory-level performance on-site, offering fast, accurate results and the ease-of-use analysts demand in the field. With a measurement range of 0 to 1000 NTU and a resolution of 0.01 NTU, the 2100P turbidimeter is ideal for regulatory monitoring, process control or field studies.

DW

ww

PW

IW

E

<del>-</del> | -

### Specifications\*

	2100P	2100P IS
<b>Measurement Method</b>	Nephelometric Ratio	
Regulatory	Meets EPA Method 180.1	Meets EN ISO 7027
Light Source	Tungsten lamp	Light-emitting diode (LED) @ 860 nm
Range		
Automatic Range Mode	0 to 1000 NTU	0 to 1000 FNU
Manual Range Selection	0 to 9.99, 0 to 99.9 and 0 to 1000 NTU	0 to 9.99, 0 to 99.9 and 0 to 1000 FNU
Accuracy	±2% of reading plus stray light	
Repeatability	±1% of reading, or 0.01 NTU, whichever is greater	±1% of reading, or 0.01 FNU, whichever is greater
Resolution	0.01 on lowest range	
Signal Averaging	Selectable on/off	
Power Requirement	4 AA alkaline batteries or optional battery eliminator	
Battery Life, Typical	300 tests with signal average mode off	
	180 tests with signal average mode on	
Operating Temperature	0 to 50°C (32 to 122°F)	
Sample Required	15 mL (0.5 oz.)	
Sample Cells	60 x 25 mm (2.36 x 1 in.) borosilicate glass with screw (	caps
Dimensions	22.2 x 9.5 x 7.9 cm (8.75 x 3.75 x 3.12 in.)	
Weight	0.5 kg (1.1 lb.); shipping weight 2.7 kg (6 lb.)	
Warranty	2 years	

\*Specifications subject to change without notice.



## **Ordering Information**

#### Hach 2100P and 2100P IS Portable Turbidimeters

Portable turbidimeters are supplied with sample cells, ampules of <0.1\*, 20, 100, and 800 NTU StablCal<sup>®</sup> Stabilized Formazin Standards, 15 mL silicone oil and sample cell oiling cloth, Instrument and Operations Manual, Quick Reference Card, four AA alkaline batteries, and carrying case.

4650000 2100P Portable Turbidimeter, USEPA compliant4474000 2100P IS Portable Turbidimeter, ISO compliant

#### StablCal® Stabilized Formazin Standards

2659400	Calibration Set, 500 mL each <0.1*, 20, 100, and 800 NTU
2659410	Calibration Set, 100 mL each <0.1*, 20, 100, and 800 NTU
2659405	Calibration Set, sealed vials <0.1*, 20, 100, and 800 NTU
2659742	<0.1* NTU, 100 mL
2659749	<0.1* NTU, 500 mL
2659753	<0.1* NTU, 1 L
2660142	20 NTU, 100 mL
2660149	20 NTU, 500 mL
2660153	20 NTU, 1 L
2660242	100 NTU, 100 mL
2660249	100 NTU, 500 mL
2660542	800 NTU, 100 mL
2660549	800 NTU, 500 mL

#### **Optional Accessories**

4397500	Degassing Kit
4397510	Sample Filtration and Degassing Kit
126936	Silicone Oil, 15 mL
4707600	Sample Cell Oiling Cloth
4607900	Battery Eliminator, 120 V
4608000	Battery Eliminator, 230 V

 $<sup>^*</sup>$ <0.1 NTU StablCal $^{@}$  Solution is used in place of low turbidity water when performing calibrations

# Other Hach Turbidimeters...

#### 2100N and 2100AN Turbidimeters

(see Lit. #2498)

Hach laboratory turbidimeters include the 2100N and 2100AN Turbidimeters, which meet USEPA Method 180.1 standards, and the 2100N IS and the 2100AN IS Turbidimeters, which meet ISO 7027 standards. All four models offer wide measurement range, ease of calibration, and the application flexibility required for demanding laboratory use.

#### **On-line Turbidimeters**

Hach also offers a complete line of advanced, on-line turbidimeters for process monitoring. Please contact Hach or your representative for more information.

Lit. No. 1619 I73 Printed in U.S.A. ©Hach Company, 2007. All rights reserved.

In the interest of improving and updating its equipment, Hach Company reserves the right to alter specifications to equipment at any time.

At Hach, it's about learning from our customers and providing the right answers. It's more than ensuring the quality of water—it's about ensuring the quality of life. When it comes to the things that touch our lives...

Keep it pure.

Make it simple.

Be right.

For current price information, technical support, and ordering assistance, contact the Hach office or distributor serving your area.

In the United States, contact:

HACH COMPANY World Headquarters P.O. Box 389

Loveland, Colorado 80539-0389

U.S.A.

Telephone: 800-227-4224 Fax: 970-669-2932 E-mail: orders@hach.com

U.S. exporters and customers in Canada, Latin America, sub-Saharan Africa, Asia, and Australia/New Zealand, contact:

HACH COMPANY World Headquarters P.O. Box 389

Loveland, Colorado 80539-0389

U.S.A.

Telephone: 970-669-3050 Fax: 970-461-3939 E-mail: intl@hach.com

In Europe, the Middle East, and Mediterranean Africa, contact:

HACH + LANGE Europe Dr. Bruno Lange GmbH & Co. KG Willstätterstraße 11 D-40549 Düsseldorf GERMANY

Tel: +49 (0) 211 5288-0 Fax: +49 (0) 211 5288-143 E-mail: info@hach-lange.de www.hach-lange.com



# TTT Environmental

Instruments and Supplies

The preferred source for instrument Rentals, Sales, Service and Supplies!

4201 B Street Anchorage, AK 99503 (907) 770-9041

www.tttenviro.com

25420 74<sup>th</sup> Avenue S. Kent, WA 98032 (253) 373-9041

info@tttenviro.com